

# Product Catalog 2019-2020

## 产品目录

# C 目录

# Contents

<b>1</b>	超低温冻存系列 (FR/FB/FT) Cryogenic Series	05-18
<b>2</b>	不锈钢耗材系列 (RTA) Stainless Steel Series	19-24
<b>3</b>	有机玻璃系列 (RTB) Acrylic Series	25-29
<b>4</b>	液体处理系列 (RTC) Liquid Handling Series	30-35
<b>5</b>	样品离心与浓缩 (RTD) Sample Centrifugation and Concentration	36-42
<b>6</b>	微生物学与细胞培养 (RTE) Microbiology & Tissue Culture	43-48
<b>7</b>	病理诊断耗材系列 (RTF) Pathology & Histology Series	49-54
<b>8</b>	样品采集系列 (RTG) Specimen Collection Series	55-57
<b>9</b>	手套系列 (RTH) Disposable Gloves Series	58-59

## 超低温冻存系列 Cryogenic Series (FR/FT/FB)

不锈钢冻存架简介 (冻存盒存储架) Ultra Low Temperature Freezer Racks Overview	06
平放抽屉型冻存架 FR-1000系列 Upright Freezer Drawer Racks (FR-1000 Series)	07
侧取抽屉型冻存架 FR-2000系列 Upright Freezer Drawer Racks (FR-2000 Series)	08
分格型冻存架 FR-3000系列 Upright Freezer Racks (FR-3000 Series)	09
单体式冻存架 FR-4000系列 Vertical Racks for Chest Freezers & Liquid Nitrogen Tank (FR-4000 Series)	10
超低温冰箱适配表 Common Freezer Matching Table	11
小型液氮罐专用吊架 FR-5000系列 Vertical Racks for Liquid Nitrogen Tank (FR-5000 Series)	12
PP冻存盒 Cryogenic Storage Boxes-PP	13
PC冻存盒 Cryogenic Storage Boxes-PC	14
纸质冻存盒 Cardboard Freezer Boxes	15
彩色翻盖纸质冻存盒 Colorful Hinged Lid Cardboard Freezer Boxes	16
冻存管 Cryogenic Tube	17
低温手套 Ultra Low Temperature Protective Glove	18

## 不锈钢耗材系列 Stainless Steel Series (RTA)

不锈钢试管架&离心管架 Stainless Steel Test Tube Racks & Centrifuge Tube Racks	20
培养皿灭菌桶 Stainless Steel Petri Dish Sterilize Container	21
不锈钢清洗篮 Stainless Steel Cleaning Basket	22
不锈钢染色架 Stainless Steel Slide Staining Racks	23
不锈钢酒精灯&涂布棒 Stainless Steel Alcohol Lamp&Stainless Steel Bacterial Cell Spreader	24

## 有机玻璃系列 Acrylic Series (RTB)

有机玻璃试管架&离心管架 Acrylic Test Tube Racks & Centrifuge Tube Racks	26
移液器架 Pipette Stand	27
移液管架 Transfer Pipette Stand	28
血凝板 Blood Clot Plate	29



## 液体处理系列 Liquid Handling Series (RTC)

吸头 Pipette Tip	31 - 32
吸头盒 Pipette Tip Box	33
巴氏吸管&试剂槽 Pasteur Pipette & Solution Basin	34
血清移液管 Serological Pipette	35

## 样品离心与浓缩 Sample Centrifugation and Concentration (RTD)

离心管 Centrifuge Tube	38 - 39
离心管盒&冰盒 Centrifuge Tube Box	40
试管 Test Tube	41
塑料试管架 Test Tube Rack	42

## 微生物学与细胞培养 Microbiology & Tissue Culture (RTE)

培养皿 Culture Dish	44
培养板 Culture Plate	45
酶标板&深孔板 Elisa Plate & Deep Well Plate	46
PCR管&PCR板 PCR Tube & PCR Plate	47
细胞刮刀&接种环 Cell Scraper & Inoculation Loop	48

## 病理诊断耗材系列 Pathology & Histology Series (RTF)

载玻片 Microscope Slides	50
盖玻片&包埋盒底模 Cover glass & Base Mold	51
载玻片盒 Slide Boxes	52
纸质邮寄夹&塑料晾片板 Cardboard Slide Mailer & Plastic Microscope Slide Tray	53
包埋盒 Embedding Cassette	54

## 样品采集系列 Specimen Collection Series (RTG)

采样杯&拭子 Urine/Stool Container & Swab	56 - 57
--	---------

## 手套系列 Disposable Gloves Series (RTH)

一次性实验手套 Disposable Gloves	59
------------------------------	----

# 1 超低温冻存系列

Cryogenic Series



### 不锈钢冻存架(冻存盒存储架) Ultra Low Temperature Freezer Racks Overview

用于冰箱冰柜内生物试样、医学临床样品的隔离存储和管理，优化存储空间。适用于Thermo Fisher、Haier、Sanyo、Eppendorf、Esco、Meiling、MVE等国内外各品牌型号冰箱、冰柜及液氮容器。

- 进口SUS304不锈钢板材制成，耐腐蚀，美观结实
- 可在零下200摄氏度环境内长期储存样本
- 每层带有独立抽屉、标签框和拉手，便于识别和拿取样本
- 四种系列：平放抽屉型、侧取抽屉型、分格型和单体式冻存架
- 适用于1英寸、2英寸以及各种非标冻存盒
- 304 grade corrosion-resistant stainless steel
- Includes pull-out handles for easy removal off the shelf
- Temperature range: stable from -196°C to 121°C
- Suitable for 1", 2" and non-standard boxes
- Match with Haier/Thermo/Sanyo/Eppendorf/Esco freezers


适用的样品盒 Cryo Box Spec	立式冰箱冻存架 Upright Freezer Racks			卧式冰柜冻存架 (兼用于液氮容器) Chest Freezer Racks
	平放抽屉型冻存架 Upright Freezer Drawer Racks	侧取抽屉型冻存架 Upright Freezer Drawer Racks	分格型冻存架 Upright Freezer Racks	单体式冻存架 Vertical Type Racks
				
1"英寸标准冻存盒 133*133*36mm	FR-1100系列 FR-1100 series	FR-2100系列 FR-2100 series	FR-3100系列 FR-3100 series	FR-4100系列 FR-4100 series
2"英寸标准冻存盒 133*133*53mm	FR-1200系列 FR-1200 series	FR-2200系列 FR-2200 series	FR-3200系列 FR-3200 series	FR-4200系列 FR-4200 series
3"英寸标准冻存盒 133*133*75mm	FR-1300系列 FR-1300 series	FR-2300系列 FR-2300 series	FR-3300系列 FR-3300 series	FR-4300系列 FR-4300 series
1.8ml*100孔塑料冻存盒 143*143*53mm	FR-1400系列 FR-1400 series	FR-2400系列 FR-2400 series	FR-3400系列 FR-3400 series	FR-4400系列 FR-4400 series
5ml*36格塑料冻存盒 146*146*65mm	FR-1500系列 FR-1500 series	FR-2500系列 FR-2500 series	FR-3500系列 FR-3500 series	FR-4500系列 FR-4500 series
72孔离心管盒 180*100*57mm	FR-1600系列 FR-1600 series	FR-2600系列 FR-2600 series	FR-3600系列 FR-3600 series	FR-4600系列 FR-4600 series
100孔离心管盒 150*150*55mm	FR-1700系列 FR-1700 series	FR-2700系列 FR-2700 series	FR-3700系列 FR-3700 series	FR-4700系列 FR-4700 series
50ml纸质离心管盒 150*150*125mm	FR-1800系列 FR-1800 series	FR-2800系列 FR-2800 series	FR-3800系列 FR-3800 series	FR-4800系列 FR-4800 series



# 超低温冻存系列

## Cryogenic Series

### 平放抽屉型冻存架 FR-1000系列 Upright Freezer Drawer Racks (FR-1000 Series)




### 技术参数 Specifications

型号 Mode	容量(H*D) Capacity	冻存盒规格 Box Spec	冻存盒尺寸(mm) Box Dimension
FR-1164	6*4	1"英寸	133*133*36
FR-1174	7*4	1"英寸	133*133*36
FR-1184	8*4	1"英寸	133*133*36
FR-1244	4*4	2"英寸	133*133*53
FR-1253	5*3	2"英寸	133*133*53
FR-1254	5*4	2"英寸	133*133*53
FR-1264	6*4	2"英寸	133*133*53
FR-1333	3*3	3"英寸	133*133*75
FR-1334	3*4	3"英寸	133*133*75
FR-1823	2*3	50ml纸质离心管盒	150*150*125
FR-1824	2*4	50ml纸质离心管盒	150*150*125

- 以上为几种常规冻存盒对应配置方案。支持定制，欢迎来电咨询
- 兼容Haier、Thermo Fisher、Sanyo、Eppendorf、Esco、Meiling等国内外各品牌型号冰箱
- Above are rack specifications for regular cryo boxes. Customized products are available
- Perfectly match with Haier, Thermo Fisher, Sanyo, Eppendorf, Esco, Meiling freezers

## 侧取抽屉型冻存架 FR-2000系列

Upright Freezer Drawer Racks (FR-2000 Series)



## 技术参数 Specifications

型号 Mode	容量(H*D) Capacity	冻存盒规格 Box Spec	冻存盒尺寸(mm) Box Dimension
FR-2164	6*4	1"英寸	133*133*36
FR-2174	7*4	1"英寸	133*133*36
FR-2184	8*4	1"英寸	133*133*36
FR-2244	4*4	2"英寸	133*133*53
FR-2253	5*3	2"英寸	133*133*53
FR-2254	5*4	2"英寸	133*133*53
FR-2264	6*4	2"英寸	133*133*53
FR-2333	3*3	3"英寸	133*133*75
FR-2334	3*4	3"英寸	133*133*75
FR-2823	2*3	50ml纸质离心管盒	150*150*125
FR-2824	2*4	50ml纸质离心管盒	150*150*125

- 以上为几种常规冻存盒对应配置方案。支持定制，欢迎来电咨询
- 兼容Haier、Thermo Fisher、Sanyo、Eppendorf、Esco、Meiling等国内外各品牌型号冰箱
- Above are rack specifications for regular cryo boxes. Customized products are available
- Perfectly match with Haier, Thermo Fisher, Sanyo, Eppendorf, Esco, Meiling freezers



# 超低温冻存系列

## Cryogenic Series

### 分格型冻存架 FR-3000系列 Upright Freezer Racks (FR-3000 Series)



#### 技术参数 Specifications

型号 Mode	容量(H*D) Capacity	冻存盒规格 Box Spec	冻存盒尺寸(mm) Box Dimension
FR-3164	6*4	1"英寸	133*133*36
FR-3174	7*4	1"英寸	133*133*36
FR-3184	8*4	1"英寸	133*133*36
FR-3244	4*4	2"英寸	133*133*53
FR-3253	5*3	2"英寸	133*133*53
FR-3254	5*4	2"英寸	133*133*53
FR-3264	6*4	2"英寸	133*133*53
FR-3333	3*3	3"英寸	133*133*75
FR-3334	3*4	3"英寸	133*133*75
FR-3823	2*3	50ml纸质离心管盒	150*150*125
FR-3824	2*4	50ml纸质离心管盒	150*150*125

- 以上为几种常规冻存盒对应配置方案。支持定制，欢迎来电咨询
- 兼容Haier、Thermo Fisher、Sanyo、Eppendorf、Esco、Meiling等国内外各品牌型号冰箱
- Above are rack specifications for regular cryo boxes. Customized products are available
- Perfectly match with Haier, Thermo Fisher, Sanyo, Eppendorf, Esco, Meiling freezers

## 单体式冻存架 FR-4000系列

Vertical Racks for Chest Freezers & Liquid Nitrogen Tank (FR-4000 Series)



### 技术参数 Specifications

型号 Mode	容量(H) Capacity	冻存盒规格 Box Spec	冻存盒尺寸(mm) Box Dimension
FR-4205	5	2" 英寸	133*133*53
FR-4206	6	2" 英寸	133*133*53
FR-4207	7	2" 英寸	133*133*53
FR-4208	8	2" 英寸	133*133*53
FR-4209	9	2" 英寸	133*133*53
FR-4210	10	2" 英寸	133*133*53
FR-4211	11	2" 英寸	133*133*53
FR-4212	12	2" 英寸	133*133*53
FR-4213	13	2" 英寸	133*133*53
FR-4309	9	3" 英寸	133*133*75
FR-4310	10	3" 英寸	133*133*75
FR-4804	4	50ml纸质离心管盒	150*150*125
FR-4805	5	50ml纸质离心管盒	150*150*125

- 以上为几种常规冻存盒对应配置方案。支持定制，欢迎来电咨询
- 兼容Haier、Thermo Fisher、Sanyo、Eppendorf、Esco、Meiling等国内外各品牌型号冰箱
- Above are rack specifications for regular cryo boxes. Customized products are available
- Perfectly match with Haier, Thermo Fisher, Sanyo, Eppendorf, Esco, Meiling freezers



# 超低温冻存系列

## Cryogenic Series

**超低温冰箱适配表**  
Common Freezer Matching Table

冰箱品牌/型号 Freezer Model	冰箱位置 Part of Freezer	2" 英寸存储方案 2" Boxes Plan				1" 英寸存储方案 1" Boxes Plan				
		冻存架 Rack Cat#	冻存架数量 Rack Plan	冻存盒数量 Box Plan	9*9盒子或者10*10盒 9*9Boxes/10*10Boxes	冻存架 Rack Cat#	冻存架数量 Rack Plan	冻存盒数量 Box Plan	10*10盒子存放 10*10Boxes	
THERMO 907	整个冰箱 Whole Freezer	FR-1254	24	480	38880/48000	FR-3174	24	672	67200	
		或/or				或/or				
		FR-2254	24	480	38880/48000	FR-1174	24	672	67200	
		或/or				或/or				
		FR-3254	24	480	38880/48000					
THERMO 705	整个冰箱 Whole Freezer	FR-1254	16	320	25920/32000	FR-3174	16	448	44800	
		或/or				或/or				
		FR-2254	16	320	25920/32000	FR-1174	16	448	44800	
		或/or				或/or				
		FR-3254	16	320	25920/32000					
SANYO MDF-U54V	冰箱上两层 Above Two Layers	FR-1264	8	192	15552/19200	FR-1184	8	256	25600	
		或/or				或/or				
			FR-3264	8	192	15552/19200	FR-3184	8	256	25600
	或/or				或/or					
	冰箱下两层 Below Two Layers	FR-3254	8	160	12960/16000	FR-1184	8	224	22400	
		或/or				或/or				
		FR-1254	8	160	12960/16000	FR-3184	8	224	22400	
SANYO MDF-U74V	冰箱上两层 Above Two Layers	FR-1264	12	288	23328/19200	FR-1184	24	768	76800	
		或/or				或/or				
			FR-3264	12	288	23328/19200	FR-3184	24	768	76800
	或/or				或/or					
	冰箱下两层 Below Two Layers	FR-3254	12	240	19440/24000	FR-3184	24	768	76800	
		或/or				或/or				
		FR-1254	12	240	19440/24000	FR-3184	24	768	76800	
HAIER DW-86L626	整个冰箱 Whole Freezer	FR-3254	20	400	32400/40000	FR-3174	20	560	56000	
		或/or				或/or				
		FR-1254	20	400	32400/40000	FR-1174	20	560	56000	
HAIER DW-86L728	整个冰箱 Whole Freezer	FR-1254	24	400	38880/48000	FR-3174	24	672	67200	
		或/or				或/or				
		FR-3254	24	400	38880/48000	FR-1174	24	672	67200	
EPPENDORF U410	整个冰箱 Whole Freezer	FR-1244	15	240	19440/24000	FR-1164	15	360	36000	
		或/or				或/or				
		FR-2244	15	240	19440/24000	FR-3164	15	360	36000	
		或/or				或/or				
		FR-3244	15	240	19440/24000					
EPPENDORF U570	整个冰箱 Whole Freezer	FR-1244	25	400	32400/40000	FR-1164	25	600	6000	
		或/or				或/or				
		FR-2244	25	400	32400/40000	FR-3164	25	600	6000	
		或/or				或/or				
		FR-2244	25	400	32400/40000					



## 小型液氮罐专用吊架 FR-5000系列

Vertical Racks for Liquid Nitrogen Tank (FR-5000 Series)

### 产品特点

- 原料：高品质SUS304不锈钢
- 耐零下200摄氏度低温、耐腐蚀
- 每层带独立小抽屉
- 适用于各种型号液氮罐及液氮箱

### FEATURES

- 304 grade corrosion-resistant stainless steel
- Temperature range: stable from -196°C to 121°C
- Independent drawer for each shelf
- Fits most standard tanks of different brands



产品编号 Item No	产品规格 Specification	冻存盒尺寸 Box Dimension (mm)	每层外尺寸 External Dim of Each Shelf (mm)
FR-5004	25格*4层 25 holes/4 shelves	71*71*45	75*75*55
FR-5005	25格*5层 25 holes/5 shelves	71*71*45	75*75*55
FR-5006	25格*6层 25 holes/6 shelves	71*71*45	75*75*55
FR-5104	25格*4层 25 holes/4 shelves	75*75*53	79*84*57
FR-5105	25格*5层 25 holes/5 shelves	75*75*53	79*84*57
FR-5106	25格*6层 25 holes/6 shelves	75*75*53	79*84*57

- FR-5004, FR-5005, FR-5006包含5\*5格不锈钢小抽屉
- FR-5104, FR-5105, FR-5106需另配5\*5格塑料液氮冻存盒(FB-2125)
- 支持来图定制

- FR-5004, FR-5005, FR-5006 include 5\*5 wells stainless steel drawers
- FR-5104, FR-5105, FR-5106 can match with 5\*5 wells PC cryo boxes (FB-2125)
- Customized products are available


# 超低温冻存系列

## Cryogenic Series

### PP冻存盒 Cryogenic Storage Boxes-PP

#### 产品特点

- 高透明聚丙烯(PP)材料制成, 可反复冻融和长时间储存样本
- 产品可耐受温度范围:  $-90^{\circ}\text{C}$ 到 $121^{\circ}\text{C}$
- 耐醇类等温和有机溶剂
- 国际标准尺寸设计, 适用于标准尺寸冻存管及冻存架
- 适用于存放0.5ml、1.0ml和2.0ml冻存管和微量管
- 多种颜色可供选择, 便于区分和查找样本



#### FEATURES

- High transparency PP material
- Temperature range: stable from  $-90^{\circ}\text{C}$  to  $121^{\circ}\text{C}$
- Chemically resistant to alcohols and mild organic solvents

产品编号 Item No	规格 Specification	隔断 Description	材料 Material	外部尺寸 Dimension(mm)	包装规格 Packing
FB-1150	2"英寸	5*10格 50 wells	PP	146*77*53	100只
FB-1181	2"英寸	9*9格 81 wells	PP	130*130*50	100只
FB-1181P	2"英寸	9*9格, 带印刷 81 wells, with printing	PP	130*130*50	100只
FB-1281	2"英寸	9*9格 81 wells	PP	133*133*53	100只
FB-1281P	2"英寸	9*9格, 带印刷 81 wells, with printing	PP	133*133*53	100只
FB-1200	2"英寸	10*10格 100 wells	PP	133*133*53	100只
FB-1300	2"英寸	10*10格 100 wells	PP	143*143*53	100只
FB-1300P	2"英寸	10*10格, 带印刷 100 wells, with printing	PP	143*143*53	100只
FB-1400	2"英寸	10*10格 100 wells	PP	141*154*55	100只
FB-1136	2.5"英寸	6*6格 36 wells	PP	145*145*126	100只
FB-1236	适用15ml离心管	4*4格 16 wells	PP	145*145*65	100只
FB-1212	适用50ml离心管	6*6格 36 wells	PP	145*145*65	100只





## PC冻存盒 Cryogenic Storage Boxes-PC

### 产品特点

- 由高强度的聚碳酸酯(PC)材料制成, 可反复冻融和长时间储存样本
- 产品可用于超低温冰箱和液氮罐中储存样本
- 产品可耐受温度范围:  $-196^{\circ}\text{C}$ 到 $121^{\circ}\text{C}$
- 耐醇类等温和有机溶剂
- 盒体的缝隙和线条设计, 利于样本的快速制冷和液体的快速排出
- 国际标准尺寸设计, 适用于标准尺寸冻存管及冻存架
- 适用于存放0.5ml、1.0ml和2.0ml冻存管和微量管
- 高透明的盒盖设计, 以及盒体上清晰的数字标识, 便于观察和查找样本
- 多种颜色可供选择, 便于区分样本
- 盒体一角的独特设计, 确保样本定位的准确性

### FEATURES

- For easy inventory management of cryogenic storage vials used in cell culture, pharmaceutical research and other biotechnology labs
- Store in mechanical freezers or liquid nitrogen
- Boxes are vented, allowing rapid cooling of contents and full drainage
- Made of tough polycarbonate for use with temperatures from  $-196^{\circ}\text{C}$  to  $121^{\circ}\text{C}$
- Chemically resistant to alcohols and mild organic solvents
- Transparent lids for PC boxes have a numbered grid for quick sample identification




产品编号 Item No	规格 Specification	隔断 Description	材料 Material	外部尺寸 Dimension(mm)	包装规格 Packing
FB-2125	2"英寸	5*5格 25 wells	PC	75*75*53	50只
FB-2181	2"英寸	9*9格 81 wells	PC	133*133*53	50只
FB-2100	2"英寸	10*10格 100 wells	PC	133*133*53	50只



# 超低温冻存系列

## Cryogenic Series

### 纸质冻存盒 Cardboard Freezer Boxes



#### 产品特点

- 由耐用型硬纸板制成，根据外层纸不同分为：普通白色，覆膜和彩色纸质冻存盒
- 耐酒精等温和有机溶剂
- 温度范围：-196°C到121°C
- 可用于液氮冷藏

#### FEATURES

- Made from durable cardboard with Superior Coatings
- Chemically resistant to alcohols and mild organic solvents
- Temperature range: stable from -196°C to 121°C
- Suitable for liquid nitrogen storage

产品编号 Item No	规格 Specification	隔断 Description	涂层种类 Coating Type	颜色 Color	外部尺寸 Dimension(mm)	包装规格 Packing
FB-3181	2"英寸	9*9格 81 wells	普通涂层 Normal coating	白色 White	133*133*50	100只
FB-3100	2"英寸	10*10格 100 wells	普通涂层 Normal coating	白色 White	133*133*50	100只
FB-3100-B	2"英寸	10*10格, 大 100 wells, big	普通涂层 Normal coating	白色 White	133*133*50	100只
FB-3281	2"英寸	9*9格 81 wells	防水涂层 Water resistant coating	白色 White	133*133*50	100只
FB-3200	2"英寸	10*10格 100 wells	防水涂层 Water resistant coating	白色 White	133*133*50	100只
FB-3381	2"英寸	9*9格 81 wells	彩色涂层 Colorful coating	彩色可选 Assorted colors	133*133*50	100只
FB-3300	2"英寸	10*10格 100 wells	彩色涂层 Colorful coating	彩色可选 Assorted colors	133*133*50	100只
FB-3481	3"英寸	9*9格 81 wells	普通涂层 Normal coating	白色 White	133*133*75	100只
FB-3400	3"英寸	10*10格 100 wells	普通涂层 Normal coating	白色 White	133*133*75	100只
FB-3581	4.5"英寸	9*9格 81 wells	普通涂层 Normal coating	白色 White	133*133*115	100只
FB-3500	4.5"英寸	10*10格 100 wells	普通涂层 Normal coating	白色 White	133*133*115	100只
FB-3136	适用于国产 5ml冻存管	6*6格 36 wells	普通涂层 Normal coating	白色 White	134*134*60	100只
FB-3336	适用于 15ml离心管	6*6格 36 wells	普通涂层 Normal coating	白色 White	136*136*115	100只
FB-3116	适用于 50ml离心管	4*4格 36 wells	普通涂层 Normal coating	白色 White	150*150*125	100只

● FB-3149为15ml离心管纸盒；FB-3116为50ml离心管纸盒


● FB-3336为15ml离心管纸盒；FB-3116为50ml离心管纸盒

● FB-3149 is 15ml centrifuge tube box; FB-3116 is 50ml centrifuge tube box

● FB-3336: 15ml centrifuge tube box; FB-3116: 50ml centrifuge tube box

## 彩色翻盖纸质冻存盒

Colorful Hinged Lid Cardboard Freezer Boxes



### 产品特点

- 覆膜防水纸质冻存盒，可反复冻融
- 盒盖带网格坐标，盒体有编号印刷，方便定位样本
- 格挡分为塑料和纸质(若放在液氮中，需将格挡换成纸质)
- 翻盖式的设计，保证了盒盖方向的唯一性，避免了搞混方向和错拿盖子的情况
- 多种颜色方便区分样品；有2英寸100格和81格，1英寸100格和81格可供选择
- 适用于放置0.5ml、1.5ml或2.0ml微型管和冻存管

### FEATURES

- The lid, bottom, and four sides feature alphanumeric identifications for fast, easy, and error-free identification of individual samples
- 5 Attractive colors (Red, Yellow, Blue, Green and White) provide preference and identification of samples
- Water resistant coating allows for improved durability
- Hinged lid prevents accidental misplace or loss

产品编号 Item No	规格 Specification	隔断 Description	隔断材料 Dividers Material	颜色 Color	外部尺寸 Dimension(mm)	包装规则 Packing
FB-4181	2"英寸	9*9格 81 wells	纸板 Cardboard	彩色可选 Assorted Colors	133*133*50	100只
FB-4100	2"英寸	10*10格 100 wells	纸板 Cardboard	彩色可选 Assorted Colors	133*133*50	100只
FB-4281	2"英寸	9*9格 81 wells	PP	彩色可选 Assorted Colors	133*133*50	100只
FB-4200	2"英寸	10*10格 100 wells	ABS	彩色可选 Assorted Colors	133*133*50	100只
FB-4600	1"英寸	10*10格 100 wells	纸板 Cardboard	彩色可选 Assorted Colors	133*133*36	100只
FB-4700	1"英寸	10*10格 100 wells	ABS	彩色可选 Assorted Colors	133*133*36	100只



# 超低温冻存系列

## Cryogenic Series

### 冻存管 Cryogenic Tube

#### 产品特点

- 由高强度、高韧性的医疗级聚丙烯制成，可反复冻融
- 最大储液量为最大刻度的80%
- 工作温度范围：-196℃到121℃，可耐受高温高压
- 盖中有密封圈设计，避免了样本渗漏，确保管内样本的安全
- 多种颜色的盖子可供选择，易于辨认区分样本

#### FEATURES

- Medical grade polypropylene, repeated freeze-thaw cycle
- Uniform, transparent tube wall, with erectable skirt
- Maximum liquid storage is 80% of the maximum scale
- Temperature resistance range: -196°C~121°C
- A variety of colors of lid, easy to identify samples



产品编号 Item No	容量 Volume	规格 Specification	材料 Material	包装规格 Packing Info
FT-05-01	0.5ml	可站立 Self Standing	PP	500PCS/bag*20
FT-10-01	1.0ml	可站立 Self Standing	PP	500PCS/bag*20
FT-15-01	1.5ml	锥形底 Conical Bottom	PP	500PCS/bag*20
FT-15-02	1.5ml	锥形底连盖 Conical Bottom	PP	500PCS/bag*20
FT-15-03	1.5ml	可站立 Self Standing	PP	500PCS/bag*20
FT-20-01	2.0ml	可站立 Self Standing	PP	500PCS/bag*20
FT-50-0101	5.0ml	模刻, 可站立 Mold Graduation, Self Standing	PP	200PCS/bag*20
FT-50-0102	5.0ml	印刷, 可站立 Printed, Self Standing	PP	200PCS/bag*20
FT-100-01	10ml	模刻, 可站立 Mold Graduation, Self Standing	PP	100PCS/bag*20



产品编号 Item No	容量 Volume	规格 Specification	材料 Material	包装规格 Packing Info
FT-10-07	1.0ml	外旋管, 可站立 External Threaded Cap, Self Standing	PP	100PCS/bag*20
FT-10-08	1.0ml	内旋管, 可站立 Internal Threaded Cap, Self Standing	PP	100PCS/bag*20
FT-15-06	1.5ml	可站立 Self Standing	PP	100PCS/bag*20
FT-20-05	2.0ml	可站立 Self Standing	PP	100PCS/bag*20
FT-20-06	2.0ml	可站立, 带印刷 Self Standing, with printing	PP	100PCS/bag*20
FT-20-07	2.0ml	外旋管, 可站立 External Threaded Cap, Self Standing	PP	100PCS/bag*20
FT-20-08	2.0ml	内旋管, 可站立 Internal Threaded Cap, Self Standing	PP	100PCS/bag*20
FT-35-07	3.5ml	外旋管, 可站立 External Threaded Cap, Self Standing	PP	100PCS/bag*20
FT-35-08	3.5ml	内旋管, 可站立 Internal Threaded Cap, Self Standing	PP	100PCS/bag*20
FT-50-07	5.0ml	外旋管, 可站立 External Threaded Cap, Self Standing	PP	100PCS/bag*20
FT-50-08	5.0ml	内旋管, 可站立 Internal Threaded Cap, Self Standing	PP	100PCS/bag*20



## 低温手套

Ultra Low Temperature Protective Glove

## 耐低温防冻手套 Cryogenic Gloves

## 产品特点

- 表面采用优质防水头层牛皮，内里采用进口防寒海绵夹层、堪培拉衬里
- 手套长度有35cm/40cm/45cm/72cm可选（长度不同价格不等）
- 适用温度范围为-10℃至-86℃
- 本产品有优异的防寒防冻保护性能



产品编号	适用范围	长度 (cm)	颜色	包装规格
FG-0001	-10℃至-86℃	35	白色	1双/包
FG-0002	-10℃至-86℃	40	白色	1双/包
FG-0003	-10℃至-86℃	45	白色	1双/包
FG-0004	-10℃至-86℃	72	白色	1双/包

## 防液氮手套 Liquid Nitrogen Protective Glove

## 产品特点

- 手套由相同的三层组成：其中两层是由绝缘材料制成的薄层，其边缘粘合在一起，因此截留了大量的绝缘空气，而无需增加额外的重量或体积；内层材料具有高保温性，使其通过毛细作用带走潮气
- 手套长度有32cm/38cm/48cm/68cm可选（长度不同价格不等）
- 适用于抵抗极度寒冷，适用温度范围为-168℃至-250℃
- 非常轻便、柔软、耐用、清洁，极其灵活，佩戴起来不会觉得沉重
- 应用于生物医学、实验室研究、工业、防液氮飞溅、冷冻食品加工LNG加注干冰研磨和其他任何极寒工作环境的防冻保护



产品编号	适用范围	长度 (cm)	颜色	包装规格
FG-1001	-168℃至-250℃	32	蓝色	1双/包
FG-1002	-168℃至-250℃	38	蓝色	1双/包
FG-1003	-168℃至-250℃	48	蓝色	1双/包
FG-1004	-168℃至-250℃	68	蓝色	1双/包



# 2 不锈钢耗材系列

Stainless Steel Series





### 不锈钢试管架&离心管架

Stainless Steel Test Tube Racks & Centrifuge Tube Racks

#### 产品特点

- 由优质不锈钢材质制成
- 耐高温高压灭菌和化学腐蚀
- 有钢板和钢丝两种类型可供选择
- 不锈钢试管架&离心管架(钢板)为圆孔设计  
不锈钢试管架&离心管架(钢丝)为方孔设计，丝网编制底
- 宽底设计，稳固不易倾倒
- 多种尺寸和孔径大小可供选择



#### FEATURES

- Made of high quality stainless steel
- Autoclavable and resistant to corrosion
- Two types to choose: steel plate and steel wire  
Stainless steel test tube racks & centrifugal tube racks (steel plate) with round hole;  
Stainless steel test tube racks & centrifugal tube racks (steel wire) with square hole and wire constructed base
- Wide base provides stability and prevents tipping
- Different dimension and size of hole optional



产品编号 Item No	品名 Description	材料 Material	孔径 Diameter(mm)	孔数 Holes	包装规格 Packing
RTA-1100	不锈钢试管架&离心管架 Stainless Steel Test Tube Racks & Centrifuge Tube Rack	钢板 Stainless Steel Plate	13/16/17/19/21/23/26	40	50 PCS
RTA-1200	不锈钢试管架&离心管架 Stainless Steel Test Tube Racks & Centrifuge Tube Rack	钢板 Stainless Steel Plate	30	6/8/10/12/24	50 PCS
RTA-2100	不锈钢试管架&离心管架 Stainless Steel Test Tube Racks & Centrifuge Tube Rack	钢丝 Stainless Steel Wire	13/16/17/19/21/23/26	40	50 PCS
RTA-2200	不锈钢试管架&离心管架 Stainless Steel Test Tube Racks & Centrifuge Tube Rack	钢丝 Stainless Steel Wire	30	6/8/10/12/24	50 PCS

注：钢丝试管架规格可定制 Customized is available

# 不锈钢耗材系列

## Stainless Steel Series

### 培养皿灭菌桶

Stainless Steel Petri Dish Sterilize Container



#### 产品特点

- 由高品质耐腐蚀不锈钢材质制成，经久耐用
- 由三部分组成：内架、桶体、盖子
- 盖与筒的上部分别有出气口，灭菌时对准，参差时密封，隔绝污染源
- 一次可以装10到20个培养皿
- 规格齐全，可定制

#### FEATURES

- High quality stainless steel material, durable
- Consists of three parts: shelf, barrel and lid
- Can hold 10pcs to 20 pcs of petri dishes at a time
- Instead of the traditional kraft paper, use simply and quickly
- Specifications can be customized

产品编号 Item No	品名 Description	规格 Specification	尺寸 Dimension	包装规格 Packing
RTA-3001	60mm平皿灭菌桶 60mm Stainless Steel Petri Dish Sterilize Container	10 PCS/Container	Φ75*180	30 PCS
RTA-3002	75mm平皿灭菌桶 75mm Stainless Steel Petri Dish Sterilize Container	10 PCS/Container	Φ90*210	30 PCS
RTA-3003	90mm平皿灭菌桶 90mm Stainless Steel Petri Dish Sterilize Container	10 PCS/Container	Φ108*230	30 PCS
RTA-3004	90mm平皿灭菌桶 90mm Stainless Steel Petri Dish Sterilize Container	12 PCS/Container	Φ108*250	30 PCS

注：桶高可按照客户需求定制 Customized is available

产品编号 Item No	品名 Description	尺寸 Dimension	包装规格 Packing
RTA-4001	移液管灭菌桶 Stainless Steel Pipette Sterilize Container	Φ63*230	50 PCS
RTA-4002	移液管灭菌桶 Stainless Steel Pipette Sterilize Container	Φ63*350	50 PCS
RTA-4003	移液管灭菌桶 Stainless Steel Pipette Sterilize Container	Φ75*350	50 PCS
RTA-4004	移液管灭菌桶 Stainless Steel Pipette Sterilize Container	Φ90*350	50 PCS
RTA-4005	移液管灭菌桶 Stainless Steel Pipette Sterilize Container	Φ108*350	50 PCS

注：桶高可按照客户需求定制 Customized is available




### 不锈钢清洗篮

Stainless Steel Cleaning Basket

广泛使用在实验室各类器皿、器械零部件、手术室器械等清洗及高压蒸汽灭菌。  
Widely used in cleaning and autoclaving of labware and medical apparatus.

#### 产品特点

- 由优质SUS304不锈钢材料制成
- 焊接牢固，无突出焊点，无毛边，无脱落，使用安全
- 通体光亮，不生锈，美观卫生
- 无毒，环保，耐高温，耐酸碱，耐腐蚀
- 网格化设计利于水或蒸汽的穿透，保证清洗和灭菌效果



#### FEATURES

- Made of high quality stainless steel
- Firm welding, no burrs, no loss, and use safely
- Perfect light, and do not rust
- Non-toxic, environmental protection, autoclavable and resistant to corrosion
- Grid design is conducive to water or steam penetration, which ensures the effect of cleaning and sterilization



产品编号 Item No	品名 Description	规格 Specification	材料 Material	包装规格 Packing
RTA-5000	不锈钢清洗篮 Stainless Steel Cleaning Basket	定制 Customized	SUS304	1 PCS

注：尺寸可按照客户需求定制 Customized is available

# 不锈钢耗材系列

Stainless Steel Series

## 不锈钢染色架 Stainless Steel Slide Staining Rack



### 产品特点

- 由优质不锈钢材料制成
- 耐腐蚀，不易生锈
- 有提手，使用更加便利
- 容易清洗
- 经久耐用
- 按照规格不同，分别可放置26片、30片和60片玻片

### FEATURES

- High quality stainless steel material, durable
- Resistant to corrosion
- Handles allow for easy carrying
- 3 Types:26 pcs,30 pcs and 60 pcs

产品编号 Item No	品名 Description	规格 Specification	包装规格 Packing
RTA-6001	不锈钢染色架 Stainless Steel Slide Staining Rack	26pcs	50 PCS
RTA-6002	不锈钢染色架 Stainless Steel Slide Staining Rack	30pcs	50 PCS
RTA-6003	不锈钢染色架 Stainless Steel Slide Staining Rack	60pcs	50 PCS



### 不锈钢酒精灯

Stainless Steel Alcohol Lamp

广泛用于实验室、工厂、医疗、科研行业等，是以酒精为燃料的加热工具。



RTA-9001

#### 产品特点

- 由优质不锈钢材料制成
- 焊接牢固，更具安全性
- 外表抛光处理，耐腐蚀，便于清洗、消毒
- 耐热性好且性能稳定，使用寿命长
- 外形美观、精巧，操作简单，便于携带
- 酒精灯整体高度低、底部宽大，可防止酒精灯倾倒溢出，大大减少事故发生


#### FEATURES

- Made of high quality stainless steel
- Firm welding, no burrs, no loss
- Perfect light, resistant to corrosion
- Beautiful design, easy to operate and carry
- Low overall height and wide bottom greatly reduce the accident

### 涂布棒

Stainless Steel Bacterial Cell Spreader

涂布棒是实验室常用品，用于培养细菌时涂板。



RTA-9002

#### 产品特点

- 由优质不锈钢材料制成
- 可高压灭菌
- 表面平滑，操作简单
- 坚固耐用，不像玻璃涂布棒那样易碎，使用寿命更长
- 便于携带

#### FEATURES

- Made of high quality stainless steel
- Autoclavable and resistant to corrosion
- Smooth surface, and simple operation
- Strong and durable
- Easy to carry

产品编号 Item No	品名 Description	规格 Specification	容量/尺寸 Vol/Dim	包装规格 Packing
RTA-9001	不锈钢酒精灯 Stainless Steel Alcohol Lamp	200ml	200ml	50PCS
RTA-9002	涂布棒 Stainless Steel Bacterial Cell Spreader	三角 trilateral	35mm*200mm (柄长) 35mm*200mm (length)	100PCS

# 3 有机玻璃系列

Acrylic Series






### 有机玻璃试管架&离心管架

Acrylic Test Tube Racks & Centrifuge Tube Racks

分别设计用于盛放试管和离心管



#### 产品特点

- 由优质有机玻璃（PMMA）材质制成
- 半透明，耐化学腐蚀
- 极佳的耐气候型，可在室外使用
- 保养和清洗方便
- 多种尺寸可供选择

#### FEATURES

- Made of high quality PMMA
- Translucent, resistance to most chemicals
- Excellent weather resistance, can be used outdoor
- Easy to maintain and clean
- Available with various size


产品编号 Item No	品名 Description	规格 Specification	材料 Material
RTB-1001	离心管架Centrifuge Tube Rack	0.5ml*40孔/holes	PMMA
RTB-1002	离心管架Centrifuge Tube Rack	1.5ml*24孔/holes	PMMA
RTB-1003	离心管架Centrifuge Tube Rack	5ml*24孔/holes	PMMA
RTB-1004	离心管架Centrifuge Tube Rack	7ml*24孔/holes	PMMA
RTB-1005	离心管架Centrifuge Tube Rack	10ml*18孔/holes	PMMA
RTB-1006	离心管架Centrifuge Tube Rack	15ml*18孔/holes	PMMA
RTB-1007	离心管架Centrifuge Tube Rack	50ml*8孔/holes	PMMA
RTB-1008	离心管架Centrifuge Tube Rack	50ml*10孔/holes	PMMA
RTB-1009	离心管架Centrifuge Tube Rack	100ml*6孔/holes	PMMA
RTB-1010	离心管架Centrifuge Tube Rack	100ml*8孔/holes	PMMA
RTB-1011	试管架Test Tube Rack	Φ13*40孔/holes	PMMA
RTB-1012	试管架Test Tube Rack	Φ16*40孔/holes	PMMA
RTB-1013	试管架Test Tube Rack	Φ19*40孔/holes	PMMA
RTB-1014	试管架Test Tube Rack	Φ21*40孔/holes	PMMA
RTB-1015	试管架Test Tube Rack	Φ23*40孔/holes	PMMA
RTB-1016	试管架Test Tube Rack	Φ26*40孔/holes	PMMA

# 有机玻璃系列

## Acrylic Series

### 移液器架 Pipette Stand

适用于放置移液器，使移液器在取用时不接触底端，从而避免潜在污染的可能性。



#### 产品特点

- 由优质有机玻璃（PMMA）材质制成
- 结实稳定，使用方便
- 移液器架有Z型和圆盘两种型号可供选择  
Z型可放置5支或6支移液器，有白色（常规）、蓝色、红色三种颜色  
圆盘可放置6支移液器，颜色为白色
- 适用于大龙、百得、艾本德等各种国产和进口移液器

#### FEATURES

- Made of high quality PMMA
- Strong stability, and easy to use
- Two models available: "Z" type and disc type  
"Z" type can place five or six pipettes;  
Disc type can place six pipettes
- Suitable for a variety of domestic and imported pipettes

产品编号 Item No	品名 Description	规格 Specification	多色可选 Assorted Colors	材料 Material	包装规格 Packing
RTB-2001-5	移液器架 Pipette Stand	Z型（5 & 6支装） Z type (5/6pieces)	白色/蓝色/红色 White/blue/red	有机玻璃	20只
RTB-2001-6	移液器架 Pipette Stand	Z型（6支装） Z Type (6 pieces)	白色 White	有机玻璃	20只
RTB-2002	移液器架 Pipette Stand	圆盘（6支装） Disc Type (6 pieces)	白色 White	有机玻璃	20只
	塑料移液器架	"L型"底	彩色	聚苯乙烯	20只
	塑料移液器架	"三角形"底	彩色	聚苯乙烯	20只

注：有机玻璃移液器架规格可定制 Customized is available



### 移液管架 Transfer Pipette Stand

适用于放置移液管，可以使移液管摆放整齐，取放方便，移液管之间不会发生碰撞，移液管内残留液体不会流出架外污染环境。

#### 产品特点

- 由优质有机玻璃（PMMA）材质制成
- 结实稳定，使用方便
- 清洁、美观，节省空间
- 移液管架有梯型和圆盘两种型号可供选择  
梯型移液管架可放置12支移液管  
圆盘移液管架可拆卸，便于运输和存储，可放置24支移液管

#### FEATURES

- Made of high quality PMMA
- Strong stability, and easy to use
- Clean and beautiful, save the space
- Two models available: ladder type and disc type  
Ladder type stand can place 12pcs of transfer pipette  
Disc type stand is removable, easy to transport and storage, can place 24pcs of transfer pipette




产品编号 Item No	品名 Description	规格 Specification	颜色 Colors	材料 Material	包装规格 Packing
RTB-3001	移液管架 Transfer Pipette Stand	梯型（12支装） Ladder Type (12 pieces)	White	PMMA	20只
RTB-3002	移液管架 Transfer Pipette Stand	圆型（24支装） Circular Type (24 pieces)	White	PMMA	20只
RTB-3003	移液管架 Transfer Pipette Stand	梯型（14支装） Ladder Type (14 pieces)	透明 Transparent	PMMA	20只

# 有机玻璃系列

## Acrylic Series

### 血凝板 Blood Clot Plate

适用于血凝和血凝抑制实验



#### 产品特点

- 由优质有机玻璃（PMMA）材质制成
- 可多次重复使用
- 孔壁光滑透亮，便于观察
- 按孔数不同，可分为96孔和24孔
- 按底部不同，可分为V型和U型底
- 另有一次性血凝板供应，避免清洗，使用方便

#### FEATURES

- Made of high quality PMMA
- Can be used repeatedly
- Smooth and bright hole wall, convenient to observe
- 96-well and 24-well plates optional
- V-shape bottom and U-shape bottom plates optional
- Disposable blood clot plates available, sterilized and easy to use

产品编号 Item No	品名 Description	规格 Specification	材料 Material	包装规格 Packing
RTB-4001	血凝板 Blood Clot Plate	96孔 V型90°/110° 96 well V-shape 90°/110°	PMMA	2PCS/bag*25
RTB-4002	血凝板 Blood Clot Plate	96孔U型 96 well U-shape	PMMA	2PCS/bag*25
RTB-4003	血凝板 Blood Clot Plate	24孔U型 24 well U-shape	PMMA	2PCS/bag*25
RTB-4004	血凝板 Blood Clot Plate	一次性血凝板 V型 Disposable V-shape	PVC	10PCS/bag*200
RTB-4005	血凝板 Blood Clot Plate	一次性血凝板 U型 Disposable U-shape	PVC	10PCS/bag*200

### 血清瓶&试剂瓶架 Serum Bottle & Reagent Bottle Rack

适用于放置各种规格的血清瓶和试剂瓶，取放方便，避免污染。

产品编号 Item No	品名 Description	规格 Specification	材料 Material	包装规格 Packing
RTB-5000	血清瓶&试剂瓶架 Serum Bottle & Reagent Bottle Rack	100ml/250ml/500ml/1000ml	PMMA	20PCS


# 4 液体处理系列

Liquid Handling Series





### 吸头 Pipette Tip



#### 产品特点


- 由优质医疗级聚丙烯(PP)材料制成，可耐受高温高压
- 10 $\mu$ l、200 $\mu$ l、300 $\mu$ l、1000 $\mu$ l、5ml、10ml五种容量可供选择
- 产品为通用型设计，可与大多数品牌和型号的移液器配套
- 独特的工艺设计减少了液体挂壁情况的发生
- 产品无DNase，无RNase，无致热源
- 不同的颜色便于区分
- 移液处理快速准确

#### FEATURES

- Medical grade PP material, autoclavable
- Sizes include 10 $\mu$ l, 200 $\mu$ l, 300 $\mu$ l, 1000 $\mu$ l, 5ml, 10ml
- Fits most brands and models of pipettors
- Rnase & Dnase Free, Non-pyrogenic
- Various color optional easy to identify
- Fast and accurate pipetting processing

# 液体处理系列

## Liquid Handling Series



产品编号 Item No	品名	Description	容量 Volume	材料 Material	包装规格 Packing
RTC-0010-011	配吉尔森移液器	Tip for Gilson	10 $\mu$ l	PP	1000PCS/bag*70
RTC-0200-011	配吉尔森移液器	Tip for Gilson	200 $\mu$ l	PP	1000PCS/bag*50
RTC-0200-012	配吉尔森移液器(带刻度)	Tip for Gilson(with graduation)	200 $\mu$ l	PP	1000PCS/bag*50
RTC-1000-011	配吉尔森移液器	Tip for Gilson	1000 $\mu$ l	PP	500PCS/bag*30
RTC-1000-012	配吉尔森移液器	Tip for Gilson	1000 $\mu$ l	PP	500PCS/bag*30
RTC-0010-021	配艾本德移液器	Tip for Eppendorf	10 $\mu$ l加长	PP	1000PCS/bag*70
RTC-0200-021	配艾本德移液器	Tip for Eppendorf	200 $\mu$ l	PP	1000PCS/bag*50
RTC-0200-022	配艾本德移液器	Tip for Eppendorf	200 $\mu$ l加长	PP	1000PCS/bag*50
RTC-0300-021	配艾本德移液器	Tip for Eppendorf	300 $\mu$ l	PP	1000PCS/bag*50
RTC-1000-021	配艾本德移液器	Tip for Eppendorf	1000 $\mu$ l	PP	500PCS/bag*30
RTC-449	配求精 III 型紫色移液器	Tip for Qiujing-III	4X49	PP	1000PCS/bag*50
RTC-650	配求精 IV 型黑色移液器	Tip for Qiujing-IV	6X50	PP	1000PCS/bag*50
RTC-747	配求精 III, IV 型移液器	Tip for Qiujing-III, IV	7.4X70	PP	500PCS/bag*30
RTC-5000-01	配吉尔森移液器	Tip for Gilson	5ml	PP	300PCS/bag*14
RTC-5000-02	配吉尔森/艾本德移液器(带刻度)	Tip for Gilson/Eppendorf(with graduation)	5ml	PP	100PCS/bag*10
RTC-5000-03	配芬兰移液器(小口)	Tip for Finnipipette	5ml	PP	300PCS/bag*14
RTC-5000-04	配大龙移液器(大口)	Tip for Dalong	5ml	PP	300PCS/bag*14
RTC-10000-01	配大龙/热电移液器	Tip for Dalong/Thermo	10ml	PP	100PCS/bag*25
RTC-10000-02	配荣泰移液器	Tip for Rongtai	10ml	PP	100PCS/bag*25



产品编号 Item No	品名	Description	容量 Volume	材料 Material	包装规格 Packing
RTC-0010-0901	配吉尔森移液器带滤芯	Tip for Gilson (with filter)	10 $\mu$ l	PP	1000PCS/bag*70
RTC-0010-0902	配艾本德移液器带滤芯	Tip for Eppendorf (with filter)	10 $\mu$ l	PP	1000PCS/bag*70
RTC-0200-0901	配吉尔森移液器带滤芯	Tip for Gilson (with filter)	200 $\mu$ l	PP	1000PCS/bag*50
RTC-1000-0901	配吉尔森移液器带滤芯	Tip for Gilson (with filter)	1000 $\mu$ l	PP	500PCS/bag*30
RTC-5000-0902	配吉尔森/艾本德移液器带滤芯(带刻度)	Tip for Gilson/Eppendorf (with filter, with graduation)	5ml	PP	100PCS/bag*10



# 液体处理系列

## Liquid Handling Series


### 吸头盒 Pipette Tip Box

#### 产品特点

- 由优质聚丙烯 (PP) 材质制成，耐高温高压
- 可重复使用
- 卡扣锁便于盖紧盒子，且易于单手操作

#### FEATURES

- Molded from PP material, autoclavable
- Recyclable
- Snap-lock clasp keeps the lid securely fastened, and easy for one-hand operation



产品编号 Item No	品名 Description	容量 Volume	材料 Material	包装规格 Packing
RTCB-0010-96	吸头盒 Pipette Tip Box	10 $\mu$ l*96孔/holes	PP	120PCS
RTCB-0200-96-1	吸头盒 Pipette Tip Box	200 $\mu$ l*96孔/holes	PP	100PCS
RTCB-0200-96-2	吸头盒 Pipette Tip Box	200 $\mu$ l*96孔/holes	PP	100PCS
RTCB-1000-60	吸头盒 Pipette Tip Box	1000 $\mu$ l*60孔/holes	PP	100PCS
RTCB-1000-100-1	吸头盒 Pipette Tip Box	1000 $\mu$ l*100孔/holes	PP	60PCS
RTCB-1000-100-2	吸头盒 Pipette Tip Box	1000 $\mu$ l*100孔/holes	PP	60PCS



产品编号 Item No	品名 Description	容量 Vol	材料 Material	包装规格 Packing
RTCB-5000-28	吸头盒 Pipette Tip Box	5 $\mu$ l*28孔/holes	PP	60PCS
RTCB-10000-15	吸头盒 Pipette Tip Box	10 $\mu$ l*15孔/holes	PP	60PCS





产品编号 Item No	品名 Description	容量 Volume	材料 Material	包装规格 Packing
RTCC-0200	巴氏吸管 Pasteur Pipette	0.2ML	LDPE	1000支/包, 50包/箱
RTCC-0500	巴氏吸管 Pasteur Pipette	0.5ML	LDPE	500支/包, 80包/箱
RTCC-1000	巴氏吸管 Pasteur Pipette	1.0ML	LDPE	100支/包, 100包/箱
RTCC-3000	巴氏吸管 Pasteur Pipette	3ML	LDPE	100支/包, 100包/箱
RTCC-3001	巴氏吸管 Pasteur Pipette	3ML加长	LDPE	100支/包, 100包/箱
RTCC-5000	巴氏吸管 Pasteur Pipette	5ML	LDPE	100支/包, 50包/箱
RTCC-10000	巴氏吸管 Pasteur Pipette	10ML	LDPE	100支/包, 25包/箱



产品编号 Item No	品名 Description	容量 Volume	颜色 Color	材料 Material	包装规格 Packing
RTCD-50	试剂槽 Solution Basin	50ML	多色可选 Assorted Colors	PP	100PCS/bagx50

# 液体处理系列

## Liquid Handling Series

### 血清移液管 Serological Pipettes

#### 产品特点

- 由优质聚苯乙烯 (PS) 材质制成, 便于观测
- 减少液体挂壁, 提高取样的准确性
- 无DNA酶, 无RNA酶, 无热源
- 容量误差控制在 $\pm 2\%$ 之内

#### FEATURES

- High-grade PS material, excellent for clear observation
- Reduce liquid attachment on the pipette surface to assure accurate delivery
- Dnase & Rnase free, Non-pyrogenic
- Graduations are calibrated for accurate dispensing to within  $\pm 2\%$



产品编号 Item No	容量 Volume	长度 Length	灭菌 Sterile	刻度 Graduations	色环 Color Ring	包装规格 Packing
RTCE-0010	1ml	278mm	灭菌	1/100	黄色	独立包装, 1000支/盒, 4盒/箱
RTCE-0020	2ml	278mm	灭菌	1/100	绿色	独立包装, 800支/盒, 4盒/箱
RTCE-0050	5ml	347mm	灭菌	0.1	蓝色	独立包装, 300支/盒, 6盒/箱
RTCE-0100	10ml	343mm	灭菌	0.1	红色	独立包装, 200支/盒, 6盒/箱
RTCE-0250	25ml	310mm	灭菌	0.2	紫色	独立包装, 200支/盒, 6盒/箱
RTCE-0500	50ml	365mm	灭菌	1	黑色	独立包装, 200支/盒, 4盒/箱
RTCE-0012	1ml	278mm	灭菌	1/100	黄色	1000支/盒, 4盒/箱
RTCE-0022	2ml	278mm	灭菌	1/100	绿色	800支/盒, 4盒/箱
RTCE-0052	5ml	347mm	灭菌	0.1	蓝色	500支/盒, 6盒/箱
RTCE-0102	10ml	343mm	灭菌	0.1	红色	350支/盒, 6盒/箱
RTCE-0252	25ml	310mm	灭菌	0.2	紫色	350支/盒, 4盒/箱
RTCE-0502	50ml	365mm	灭菌	1	黑色	250支/盒, 4盒/箱

# 5

## 样品离心与浓缩

Sample Centrifugation and Concentration







# 样品离心与浓缩

## Sample Centrifugation and Concentration

### 离心管系列 Centrifuge Tube Series



#### 产品特点

- 特级PP材质，化学稳定性强
- 管体高透明，便于样本观测
- 耐高温高压灭菌(不带盖)
- 最大离心速度高达10500g(15ml)/9400g(50ml)/3000g(50ml自立式)
- 伽马射线或环氧乙烷消毒，安全无菌(同时提供非灭菌产品)
- 适用于大多数标准台式和立式离心机

#### FEATURES

- Special grade of PP material, chemical stability
- Transparent wall permits easy viewing of tube contents
- Autoclavable (without cap)
- Max.RCF rating to 10500g(15ml)/9400g(50ml)/3000g(50ml free-standing)
- Gama or EO sterilization

产品编号 Item No	品名	Description	容量 Volume	材料 Material	包装规格 Packing
RTD-100-01	离心管(圆底插口)	Centrifuge Tube(snap cap)	100ml	PP	30PCS/bagx15
RTD-100-02	离心管(圆底/螺旋盖)	Centrifuge Tube	100ml	PP	30PCS/bagx15
RTD-100-03	离心管(圆底/螺旋盖印刷)	Centrifuge Tube(printed)	100ml	PP	30PCS/bagx15
RTD-50-01	离心管(平底/螺旋盖)	Centrifuge Tube(flat bottom)	50ml	PP	50PCS/bagx15
RTD-50-02	离心管(尖底/螺旋盖)	Centrifuge Tube(conical bottom)	50ml	PP	50PCS/bagx15
RTD-50-03	离心管(圆底/螺旋盖)	Centrifuge Tube(round bottom)	50ml	PP	50PCS/bagx15
RTD-50-04	离心管(圆底/插口)	Centrifuge Tube(snap cap)	50ml	PP	50PCS/bagx15



# 样品离心与浓缩

## Sample Centrifugation and Concentration



产品编号 Item No	品名	Description	规格 Specification	材料 Material	包装规格 Packing
RTD-30-01	离心管(圆底/内旋盖)	Centrifuge Tube(internal threaded cap)	30ml	PP	50PCS/bagx15
RTD-20-01	离心管(圆底/插口)	Centrifuge Tube(snap cap)	20ml	PP	100PCS/bagx18
RTD-15-01	离心管(圆底/螺旋盖)	Centrifuge Tube(round bottom)	15ml	PP	100PCS/bagx24
RTD-15-02	离心管(尖底/螺旋盖)	Centrifuge Tube(conical bottom)	15ml	PP	100PCS/bagx24
RTD-10-01	离心管(圆底/连盖)	Centrifuge Tube(connect cap)	10ml	PP	200PCS/bagx14
RTD-10-02	离心管(圆底/螺旋盖)	Centrifuge Tube(round bottom)	10ml	PP	100PCS/bagx22
RTD-10-03	离心管(尖底/螺旋盖)	Centrifuge Tube(conical bottom)	10ml	PP	100PCS/bagx22



产品编号 Item No	品名	Description	规格 Specification	材料 Material	包装规格 Packing
RTD-7-01	离心管(圆底/连盖)	Centrifuge Tube	7ml	PP	200PCS/bagx20
RTD-5-01	离心管(圆底/连盖)	Centrifuge Tube	5ml	PP	300PCS/bagx20
RTD-5-02	离心管(圆底/连盖棕色)	Centrifuge Tube(brown color)	5ml	PP	300PCS/bagx20
RTD-2-01	离心管(圆底/连盖)	Centrifuge Tube	2ml	PP	500PCS/bagx25
RTD-1.5-01	离心管(尖底/连盖带齿)	Centrifuge Tube(with gear)	1.5ml	PP	500PCS/bagx25
RTD-1.5-02	离心管(尖底/连盖大盖)	Centrifuge Tube	1.5ml	PP	500PCS/bagx25
RTD-05-01	离心管(尖底/连盖)	Centrifuge Tube	0.5ml	PP	1000PCS/bagx28
RTD-02-01	离心管(尖底/凸盖)	Centrifuge Tube(Domed-shaped)	0.2ml	PP	1000PCS/bagx60
RTD-02-02	离心管(尖底/平盖)	Centrifuge Tube(Flat-top)	0.2ml	PP	1000PCS/bagx60

# 样品离心与浓缩

## Sample Centrifugation and Concentration


### 离心管盒&冰盒 Centrifuge Tube Box

#### 产品特点




- 由医疗级聚丙烯(PP)材质制成
- 耐酒精等温和的有机溶剂
- 工作范围温度：-196℃到121℃
- 盖子可移动，印有书写区域，方便标记

#### FEATURES

- Made from medical grade PP material
- Chemically resistant to alcohols and mild organic solvents
- Temperature range: stable from -196°C to 121°C
- Removable lid includes an inventory list writing area



产品编号 Item No	品名	Description	规格 Specification	材料 Material	包装规格 Packing
RTDB-1001	两用离心管架	Double-function Centrifuge Tube Rack	60-hole(0.5ml/1.5ml)	PP	150PCS
RTDB-1002	两用离心管架	Double-function Centrifuge Tube Rack	80-hole(0.5ml/1.5ml)	PP	100PCS
RTDB-1003	两用离心管架	Double-function Centrifuge Tube Rack	96-hole(0.5ml/1.5ml)	PP	100PCS
RTDB-1004	多用离心管架	Multi-purpose Centrifuge Tube Rack	0.5ml/1.5ml/15ml/50ml	PP	50PCS
RTDB-1005	离心管两用架	Centrifuge Tube Rack	15ml/50ml	PP	50PCS
RTDB-1006	离心管两用架	Centrifuge Tube Rack	15ml/50ml	PP	50PCS



产品编号 Item No	品名	Description	规格 Specification	材料 Material	包装规格 Packing
RTDB-1007	离心管盒	Centrifuge Tube Box	0.2ml/96-hole	PP	200PCS
RTDB-1008	离心管盒	Centrifuge Tube Box	0.5ml/72-hole	PP	80PCS
RTDB-1009	离心管盒	Centrifuge Tube Box	1.5ml/72-hole	PP	80PCS
RTDB-1010	离心管盒	Centrifuge Tube Box	1.5ml/100-hole	PP	80PCS
RTDB-1011	离心管盒	Centrifuge Tube Box	5ml/32-hole	PP	80PCS
RTDB-1012	离心管冰盒	Ice Box	1.5ml/50-hole	PP	80PCS
RTDB-1013	离心管冰盒	Ice Box	1.5ml/96-hole	PP	60PCS
RTDB-1014	多用离心管冰盒	Multi-purpose Ice Box	0.2ml/0.5ml/1.5ml/54-hole	PP	80PCS



# 样品离心与浓缩

## Sample Centrifugation and Concentration

### 试管 Test Tube

#### 产品特点

- 两种材料的管体可供选择：  
高透明的聚苯乙烯(PS)制成，便于样本观测；  
优质的聚丙烯(PP)制成，耐高温高压、耐化学腐蚀
- 重量轻，持久耐用，安全
- 根据底部不同可分为：圆底试管和尖底试管
- 多种规格可供选择

#### FEATURES

- Ps and PP material tube for choose
- Light weight, durable and safety
- Available with round bottom and conical bottom
- A variety of specifications to choose



产品编号 Item No	品名	Description	规格 Specification	材料 Material	包装规格 Packing
RTDC-1260-S	圆底试管	Round Bottom Test Tube	12*60mm	PS	500PCS/bag*12
RTDC-1260-P	圆底试管	Round Bottom Test Tube	12*60mm	PP	500PCS/bag*12
RTDC-1275-S	圆底试管	Round Bottom Test Tube	12*75mm	PS	500PCS/bag*10
RTDC-1275-P	圆底试管	Round Bottom Test Tube	12*75mm	PP	500PCS/bag*10
RTDC-1375-S	圆底试管	Round Bottom Test Tube	13*75mm	PS	500PCS/bag*10
RTDC-1375-P	圆底试管	Round Bottom Test Tube	13*75mm	PP	500PCS/bag*10
RTDC-1310-S	圆底试管	Round Bottom Test Tube	13*100mm	PS	400PCS/bag*10
RTDC-1310-P	圆底试管	Round Bottom Test Tube	13*100mm	PP	400PCS/bag*10
RTDC-1510-S	圆底试管	Round Bottom Test Tube	15*100mm	PS	300PCS/bag*10
RTDC-1510-P	圆底试管	Round Bottom Test Tube	15*100mm	PP	300PCS/bag*10
RTDC-1610-S	圆底试管	Round Bottom Test Tube	16*100mm	PS	250PCS/bag*10
RTDC-1610-P	圆底试管	Round Bottom Test Tube	16*100mm	PP	250PCS/bag*10
RTDC-2212	圆底试管	Round Bottom Test Tube	22*120mm	PS	100PCS/bag*12
RTDC-0005	圆底试管	Round Bottom Test Tube	5ml	PS	50PCS/bag*80
RTDC-0010	圆底试管	Round Bottom Test Tube	10ml	PS	1000PCS/bag*4
RTDC-0012-1	直口尿沉渣管	Urine Test Tube (Direct Hole)	12ml	PS	100PCS/bag*24
RTDC-0012-2	反口尿沉渣管	Urine Test Tube (Bugle Hole)	12ml	PS	100PCS/bag*24
RTDC-0012-3	螺口尿沉渣管	Urine Test Tube (with lid)	12ml	PS	100PCS/bag*24
RTDC-S-12	试管塞(各种颜色)	Tube Stopper(Various Colors)	Φ12	LDPE	1000PCS/bag*30
RTDC-S-13	试管塞(各种颜色)	Tube Stopper(Various Colors)	Φ13	LDPE	1000PCS/bag*30
RTDC-S-15	试管塞(各种颜色)	Tube Stopper(Various Colors)	Φ15	LDPE	1000PCS/bag*15
RTDC-S-16	试管塞(各种颜色)	Tube Stopper(Various Colors)	Φ16	LDPE	1000PCS/bag*15



# 样品离心与浓缩

## Sample Centrifugation and Concentration


### 塑料试管架 Test Tube Rack

#### 产品特点

- 由优质聚丙烯(PP)材质或ABS塑料制成
- 耐酒精等温和的有机溶剂
- 多功能试管架上的硅胶垫可以固定多种规格的试管

#### FEATURES

- Made from high quality PP and ABS material
- Chemically resistant to alcohols and mild organic solvents
- Silicon rubber pad lay on the rack ensure tubes of various diameters be fixed in the rack



产品编号 Item No	品名	Description	规格 Specification	材料 Material	包装规格 Packing
RTDD-1001	立柱式试管架	Pillar Stand Test Tube Rack	5*10/Φ13	PP	50PCS
RTDD-1002	立柱式试管架	Pillar Stand Test Tube Rack	5*16/Φ15	PP	50PCS
RTDD-1003	试管架	Test Tube Rack	5*10/Φ13	PP	100PCS
RTDD-1004	试管架	Test Tube Rack	5*10/Φ16	PP	100PCS
RTDD-1005	试管架	Test Tube Rack	4*10/Φ18	PP	100PCS



产品编号 Item No	品名	Description	规格 Specification	材料 Material	包装规格 Packing
RTDD-1006	试管架	Test Tube Rack	21-holes	PP	50PCS
RTDD-1007	试管架	Test Tube Rack	40-holes	PP	50PCS
RTDD-1008	试管架	Test Tube Rack	60-holes	PP	50PCS
RTDD-1009	试管架	Test Tube Rack	90-holes	PP	100PCS
RTDD-1010	多功能试管架	Multi-Function Test Tube Rack	50-holes 10-18mm	PP	100PCS

# 6 微生物学与细胞培养

Microbiology & Tissue Culture





### 培养皿 Culture Dish

#### 产品特点

- 由医疗级聚苯乙烯(PS)材料制成，产品透明度高，便于实验观察
- 产品壁厚均匀，表面光滑且不易损坏
- 上皿有三个透气孔设计，便于内外气体交换
- 经环氧乙烷灭菌

#### FEATURES

- Made of high-grade medical PS, excellent optical clarity
- Uniform and smooth surface, reliable and difficult to damage
- Three vents molded on lid, allowing for gas exchange
- EO sterilization



产品编号 Item No	品名 Description	规格 Specification	材料 Material	包装规格 Packing
RTE-35	培养皿 Petri Dish	35*15mm	PS	10只/包, 200包/箱
RTE-60	培养皿 Petri Dish	60*15mm	PS	10只/包, 200包/箱
RTE-70	培养皿 Petri Dish	70*15mm	PS	10只/包, 80包/箱
RTE-90	培养皿 Petri Dish	90*15mm	PS	10只/包, 50包/箱
RTE-90-2	分格培养皿 Petri Dish	90*15mm	PS	10只/包, 50包/箱
RTE-90-3	三格培养皿 Petri Dish	90*15mm	PS	10只/包, 50包/箱
RTE-90-4	四格培养皿 Petri Dish	90*15mm	PS	10只/包, 50包/箱
RTE-100	培养皿 Petri Dish	100*15mm	PS	10只/包, 50包/箱
RTE-120	培养皿 Petri Dish	120*15mm	PS	10只/包, 32包/箱
RTE-150	培养皿 Petri Dish	150*15mm	PS	10只/包, 20包/箱
RTE-100-2	培养皿 Petri Dish	方形, 100*100mm	PS	10只/包, 50包/箱
RTE-130	培养皿 Petri Dish	方形, 130*130mm	PS	10只/包, 30包/箱
RTE-250	培养皿 Petri Dish	方形, 250*250mm	PS	4只/包, 10包/箱



产品编号 Item No	品名 Description	规格 Specification	材料 Material	包装规格 Packing
RTE-100	培养皿 Petri Dish	100*150mm	PS	10PCS/bag*50
RTE-120	培养皿 Petri Dish	120*15mm	PS	10PCS/bag*30
RTE-150	培养皿 Petri Dish	150*15mm	PS	10PCS/bag*20
RTE-100-2	培养皿 Petri Dish	100*100mm	PS	10PCS/bag*50
RTE-130	培养皿 Petri Dish	130*130mm	PS	10PCS/bag*15
RTE-250	培养皿 Petri Dish	250*250mm	PS	1PCS/bag*40



# 微生物学与细胞培养

## Microbiology & Tissue Culture

### 培养板 Culture Plate



#### 产品特点

- 由医疗级聚苯乙烯(PS)材料制成
- 透明度高, 兼容于常规自动化系统
- 防滑底边, 方便取放
- 数字和字母编码, 易于识别
- 可堆叠

#### FEATURES

- Made of high-grade medical PS
- High transparency, compatible with common automated systems
- Anti-slip area on base sides is easy to grip
- Alpha-Numeric codes for easy identification
- Stackable



产品编号 Item No	品名 Description	规格 Specification	材料 Material	包装规格 Packing
RTEB-96-1	培养板U型底 Culture Plate(U bottom)	96孔/holes	PS	200只
RTEB-96-2	培养板V型底 Culture Plate (V bottom)	96孔/holes	PS	200只
RTEB-96-3	培养板 平底 Culture Plate (Flat bottom)	96孔/holes	PS	200只
RTEB-48	培养板 Culture Plate	48孔/holes	PS	200只
RTEB-24	培养板 Culture Plate	24孔/holes	PS	200只
RTEB-12	培养板 Culture Plate	12孔/holes	PS	200只
RTEB-06	培养板 Culture Plate	6孔/holes	PS	200只
RTEB-384-1	培养板 Culture Plate	384孔, 高款	PS	200只
RTEB-384-2	培养板 Culture Plate	384孔, 矮款	PS	200只

### 酶标板 Elisa Plate



#### 产品特点


- 由优质聚苯乙烯(PS)材质制成
- 孔数为96孔
- 平底酶标板, 适合在酶标仪检测
- 可拆式的设计, 可分为板条和外框, 板条有8连和12连可选择

#### FEATURES

- Made of high-grade medical PS
- 96-well plate
- Flat bottom, suitable for ELIASA
- Detachable, and can be divided into strips (8-stripe or 12-stripe) and the outer frame

产品编号 Item No	品名 Description	规格 Specification	材料 Material	包装规格 Packing
RTEC-9608	平底酶标板 Elisa Plate	96孔(8连可拆) 96 well(detachable 8 strips)	PS	200PCS
RTEC-9612	平底酶标板 Elisa Plate	96孔(12连可拆) 96 well(detachable 12 strips)	PS	200PCS

### 深孔板 Deep-well Plate



#### 产品特点

- 由优质聚丙烯(PP)材料制成
- 1ml和2ml的设计可以提供最适宜的容量选择
- PP材质可耐受121°C高温高压灭菌处理


#### FEATURES

- Made from high quality polypropylene (PP)
- Can provide the most suitable capacity: 1ml and 2ml
- Autoclavable to 121°C

产品编号 Item No	品名 Description	容量 Volume	材料 Material	包装规格 Packing
RTED-96-1	深孔板 Deep Well Plate	1.2ml*96well	PP	100PCS
RTED-96-2	深孔板 Deep Well Plate	1.6ml*96well	PP	100PCS
RTED-96-2	深孔板 Deep Well Plate	2.2ml*96well	PP	100PCS
RTED-96-8	深孔板盖子 Cover	—	PP	100PCS
RTED-96-9	深孔板硅胶盖 Silicon Sheet	96 well	硅胶 Silicon	10PCS/bag*20



### PCR管 PCR Tube



#### 产品特点


- 由优质医疗级聚丙烯（PP）材料制成
- 工作温度范围：-20°C到100°C
- 适用于常规品牌和型号的PCR仪
- 模注的横向纵向坐标标识，便于样本识别

#### FEATURES

- Materials: Virgin polypropylene
- Temperature range: stable from -20°C to 100°C
- Fits conventional brands and models of PCR instrument
- Horizontal and vertical coordinates to facilitate the identification of samples

产品编号 Item No	品名 Description	容量 Volume	材料 Material	包装规格 Packing
RTEE-101	薄壁PCR管 PCR Tube	0.2ml	PP	1000PCS/bag*15
RTEE-108	八联PCR管 PCR Tube	0.2ml/8 strips	PP	200PCS/bag*15
RTEE-112	十二联PCR管 PCR Tube	0.2ml/12 strips	PP	200PCS/bag*15


### PCR板 PCR Plate




产品编号 Item No	品名 Description	容量 Volume	材料 Material	包装规格 Packing
RTEE-201	PCR板 PCR Plate	96well	PP	20PCS/bag*10
RTEE-202	PCR板 PCR Plate	96孔凸 96well/ bulge	PP	20PCS/bag*10
RTEE-203	PCR板 PCR Plate	96孔带边 96well/ skirted	PP	20PCS/bag*10
RTEE-204	PCR板 PCR Plate	384well	PP	20PCS/bag*10



### 细胞刮刀 & 接种环 Cell Scraper & Inoculation Loop



产品编号 Item No	品名 Description	刀柄长度 cm Length of handle	刀片长度 cm Length of blade	材料 Material	包装规格 Packing
RTEF-901	细胞刮刀	25	1.8	ABS+TPE	1PCS/bag*200
RTEF-902	细胞铲	18	1.8	PE	1PCS/bag*200



产品编号 Item No	描述 Description	容量 Volume	消毒 Sterile	材料 Material	包装规格 Packing
RTEG-911	接种环 Inoculation Loop	1 $\mu$ l	YES	PS	20PCS*50peel pack*10 boxes
RTEG-912	接种环 Inoculation Loop	10 $\mu$ l	YES	PS	20PCS*50peel pack*10 boxes
RTEG-913	接种环 Inoculation Loop	1 $\mu$ l+10 $\mu$ l	YES	ABS	10PCS*100peel pack*10 boxes

# 7 病理诊断耗材系列

Pathology & Histology Series



### 载玻片 Microscope Slides



#### 产品特点

- 由优质玻璃原材料制成
- 玻片平面整洁光滑，确保可以更清晰的观察样本
- 玻片边沿经过打磨处理，四角经过倒角处理，使用操作更安全
- 非磨砂面玻片可用表面面积大；磨砂面玻片便于对样本进行标记和管理

#### FEATURES

- Made from super grade glass and super white grade material
- Clean and smooth slide, ensure to observe samples clearly
- Ground Edges, beveled edges provide utmost safety during manual operations
- Plain microscope slides supplied with large available surface area; Frosted microscope slides supplied with frosted areas on the slide end for marking

产品编号 Item No	品名 Description	规格 Specification	包装规格 Packing
RTF-7101	载玻片7101(免清洗) Microscope Slides 7101(Plain)	90°直角抛光边 Corners 90°, Ground Edges	50PCS/box*50
RTF-7105	载玻片7105单头单面磨砂 Microscope Slides(Frosted one end and one side)	45°直角抛光边 Corners 45°, Ground Edges	50PCS/box*50
RTF-7107	载玻片7107单头双面磨砂 Microscope Slides(Frosted one end and both sides)	45°直角抛光边 Corners 45°, Ground Edges	50PCS/box*50



# 病理诊断耗材系列

## Pathology & Histology Series


### 盖玻片 Cover Glass

#### 产品特点


- 由优质玻璃原材料制成
- 极佳的光学特性，稳定的分子结构
- 平整度和尺寸均一性好
- 可以根据需要提供特殊规格

#### FEATURES

- Super grade glass and super white grade material
- Ideal optical properties, stable molecular structure
- Uniform flatness and uniform sizes
- Special dimensions available upon request




产品编号 Item No	品名 Description	规格 Specification	包装规格 Packing
RTFB-1818	盖玻片 Cover Glass(Vacuum-package)	18*18mm	100PCS/box*500
RTFB-2020	盖玻片 Cover Glass(Vacuum-package)	20*20mm	100PCS/box*500
RTFB-2222	盖玻片 Cover Glass(Vacuum-package)	22*22mm	100PCS/box*500
RTFB-2424	盖玻片 Cover Glass(Vacuum-package)	24*24mm	100PCS/box*500
RTFB-2432	盖玻片 Cover Glass(Vacuum-package)	24*32mm	100PCS/box*200
RTFB-2450	盖玻片 Cover Glass(Vacuum-package)	24*50mm	100PCS/box*200
RTFB-2460	盖玻片 Cover Glass(Vacuum-package)	24*60mm	100PCS/box*200



产品编号 Item No	品名 Description	规格 Specification	包装规格 Packing
RTFC-0001	包埋盒底模 Base Mold for embedding cassette	7*7*7	50pcs/bag*10
RTFC-0002	包埋盒底模 Base Mold for embedding cassette	15*15*7	50pcs/bag*10
RTFC-0003	包埋盒底模 Base Mold for embedding cassette	24*24*7	50pcs/bag*10
RTFC-0004	包埋盒底模 Base Mold for embedding cassette	30*24*7	50pcs/bag*10
RTFC-0005	包埋盒底模 Base Mold for embedding cassette	37*24*7	50pcs/bag*10

### 载玻片盒 Microscope Slide Boxes



#### 产品特点

- 塑料载玻片盒和塑料晾片板由ABS材质制成，塑料邮寄夹由优质聚丙烯材质制成，结实美观，有效防虫防潮
- 可存放25mmx75mm的标准级载玻片
- 卓越的设计，防止玻片损坏和灰尘污染玻片
- 在玻片盒内部每个存放玻片位置都有相应的编号，便于查找和管理样本
- 100片的载玻片盒有金属扣设计，确保了盒体密封，有效地防止灰尘进入，可以长期保存和保护内部标本
- 多种颜色可供选择，可以更好地进行标本的颜色管理

#### FEATURES

- Plastic microscope slide boxes and plastic dry plates are made by high quality ABS material, and plastic slide mailers are made by high quality PP material. Ideal optical properties, stable molecular structure
- Accommodates standard microscope slides in dimensions of 25.0mmx75.0mm
- Excellent design, prevent microscope slide damage and pollution
- An index card mounted inside the lids allows for sample specimen identification
- Available in 100-place with metal buckle design, which ensure the sealing and effectively prevent the dust, preserve and
- Different colors to choose, can be better to manage the sample

产品编号 Item No	品名 Description	规格 Specification	材料 Material	包装规格 Packing
RTFD-1001	塑料邮寄夹 Slide Mailer	1PCS	PP	200PCS/bag*10
RTFD-1002	塑料邮寄夹 Slide Mailer	2PCS	PP	100PCS/bag*10
RTFD-1003	塑料邮寄夹 Slide Mailer	3PCS	PP	100PCS/bag*10
RTFD-1005	塑料邮寄夹 Slide Mailer	5PCS	PP	100PCS/bag*10
RTFD-1010	塑料载玻片盒 Slide Storage Box	10PCS	ABS	200PCS
RTFD-1012	塑料载玻片盒 Slide Storage Box	12PCS	ABS	200PCS
RTFD-1025	塑料载玻片盒 Slide Storage Box	25PCS	ABS	150PCS
RTFD-1050	塑料载玻片盒 Slide Storage Box	50PCS	ABS	150PCS
RTFD-1100-1	塑料载玻片盒 Slide Storage Box	100PCS	ABS	60PCS
RTFD-1100-2	塑料载玻片盒 Slide Storage Box	100PCS (With Wtainless Steel Lock)	ABS	80PCS



# 病理诊断耗材系列

## Pathology & Histology Series

### 纸质邮寄夹 Cardboard Mailer



产品编号 Item No	品名 Description	规格 Specification	材料 Material	包装规格 Packing
RTFE-2001	纸质邮寄夹 Cardboard Mailer	1PCS	Cardboard	50PCS/box*20
RTFE-2002	纸质邮寄夹 Cardboard Mailer	2PCS	Cardboard	50PCS/box*10
RTFE-2003	纸质邮寄夹 Cardboard Mailer	3PCS	Cardboard	50PCS/box*10
RTFE-2004	纸质邮寄夹 Cardboard Mailer	4PCS	Cardboard	50PCS/box*10
RTFE-2010	纸质邮寄夹 Cardboard Mailer	10PCS	Cardboard	50PCS




产品编号 Item No	品名 Description	规格 Specification	材料 Material	包装规格 Packing
RTFE-2020	纸质邮寄夹 Cardboard Mailer	20PCS	Cardboard	50PCS
RTFE-2030	纸质邮寄夹 Cardboard Mailer	30PCS	Cardboard	50PCS



产品编号 Item No	品名 Description	规格 Specification	材料 Material	包装规格 Packing
RTFE-2048	纸质邮寄夹 Cardboard Mailer	48PCS	Cardboard	50PCS
RTFE-3020	塑料晾片板 Plastic Microscope Slide Tray	20PCS	ABS	100PCS



### 包埋盒 Embedding Cassette



产品编号 Item No	产品规格 Specification	规格 Specification	材料 Material	包装规格 Packing
RTFF-0001	包埋盒 Embedding Cassette	长条连盖可拆 Histology-Strip Holes	POM	5000PCS
RTFF-0002	包埋盒 Embedding Cassette	方格连盖可拆 Biopsy-Square Holes	POM	5000PCS
RTFF-0003	包埋盒 Embedding Cassette	长条可拆连盖 Histology-Strip Holes	POM	5000PCS
RTFF-0004	包埋盒 Embedding Cassette	无盖圆孔 Standard Round Holes	POM	300PCS/bag*30

8

## 样品采集系列

Specimen Collection Series






# 样品采集系列

## Specimen Collection Series

### 采样杯 Urine and Stool Containers



#### 产品特点

- 可选配勺子，匹配多种颜色和螺口盖
- 具有良好的温度和化学试剂耐受性
- 规格有5ml, 20ml, 30ml, 40ml, 60ml, 90ml, 120ml
- 可选模具刻度或标签

#### FEATURES

- Optional spoon, with assorted colors and screw cap
- Excellent temperature and chemical resistance
- Sizes range of 5ml, 20ml, 30ml, 40ml, 60ml, 90ml, 120ml
- With or without molded graduation & label

产品编号 Item No	品名 Description	容量 Volume	材料 Material	包装规格 Packing
RTG-0001	采便杯 (螺旋盖) Stool Container with screw cap	5ml	LDPE	250PCS/bag*20
RTG-0002	采便杯 Stool Container	5ml	PS	400PCS/bag*10
RTG-0003	尿杯 Urine Container	40ml	PP	40PCS/bag*20
RTG-0004	采便杯 Stool Container	40ml	PP	40PCS/bag*20



产品编号 Item No	品名 Description	容量 Volume	材料 Material	包装规格 Packing
RTG-0005	采便杯 Stool Container	20ml	PP	100PCS/bag*20
RTG-0006	尿杯 Urine Container	30ml	PS	25PCS/bag*20
RTG-0007	采便杯 Stool Container	30ml	PS	25PCS/bag*20

# 样品采集系列

## Specimen Collection Series



产品编号 Item No	品名 Description	容量 Volume	材料 Material	包装规格 Packing
RTG-0008	采便杯 Stool Container	60ml	PP	25PCS/bag*20
RTG-0009	尿杯 Urine Container	60ml	PP	25PCS/bag*20
RTG-0010	采便杯 Stool Container	90ml	PP	20PCS/bag*20
RTG-0011	尿杯 Urine Container	90ml	PP	20PCS/bag*20
RTG-0012	采便杯 Stool Container	120ml	PP	20PCS/bag*20
RTG-0013	尿杯 Urine Container	120ml	PP	20PCS/bag*20



产品编号 Item No	品名 Description	容量 Volume	材料 Material	包装规格 Packing
RTG-6001	拭子(塑料棉签棒) Swab(plastic stick)	Φ13x150mm	PP	100PCS/bag*25
RTG-6002	拭子(木质棉签棒) Swab(wooden stick)	Φ13x150mm	PP	100PCS/bag*25

9

## 手套系列

Disposable Gloves Series





# 手套系列

## Disposable Gloves Series

### 一次性实验手套 Disposable Gloves

作为医疗检测、食品加工、美容美发、家庭护理、诊断和电子厂、无尘室、洁净室、实验室等作业的手部防护工具。

#### 一次性无粉乳胶手套

- 由纯天然乳胶制成，舒适、柔软和灵活
- 五指灵活，左右手通用
- 平滑的氯化物内里表面使穿戴更方便
- 全纹理表面，有效增强抓握摩擦力
- 厚度分布均匀，增强触觉灵敏度
- 防化性能优越，耐一定的酸碱度
- 较强的抗拉强度，耐穿刺，不易破损，弹性好，透气性强
- 淡黄色

#### 一次性无粉丁腈手套


- 由丁腈橡胶材料制成，更耐穿刺和耐化学腐蚀
- 五指灵活，左右手通用
- 不含添加剂（硫磺和化学催速剂）降低IV型超敏反应
- 非乳胶配方减少人手对天然乳胶蛋白的不良反应
- 低致敏性，穿戴更安全、舒适
- 具有抗静电，耐酸碱，耐老化和耐油性
- 手指纹路设计，方便使用者抓握物体
- 卷边腕口设计，戴脱更方便
- 深蓝色

#### 芦荟手套

- 由绿色芦荟乳胶材料制成，护肤、美容、环保
- 高等级天然乳胶将人体的过敏反应降到最低，具有极强的化学抗性
- 五指灵活，左右手通用，卷边腕口
- 光滑内表面及卷边易于穿戴
- 蛋白质含量低
- 高品质天然乳胶材料，良好的拉力强度和伸长率，避免了在穿戴时撕裂
- 独特的处理工艺，无皮肤刺激、过敏现象

#### 一次性PE手套

- 由聚乙烯（PE）材料制成
- 一次性使用，防水、防油，耐酸碱，安全、无毒、环保
- 五指灵活，左右手通用
- 易穿戴不滑手，操作灵活
- 手套表面经凹凸加工，可防粘连



产品编号 Item No	类型 type	尺寸 Dimension	颜色 Color	包装规格 Packing
RTH-1001	乳胶	S/M/L	自然淡黄色	独立包装 600双/箱
RTH-2001	丁腈	S/M/L	深蓝色	独立包装 500双/箱
RTH-3001	芦荟乳胶	S/M/L	绿色	独立包装 600双/箱
RTH-4001	PE	均码	透明	100个/包*50包



**专业创造品质 诚信铸就价值**



HAORUN MEDICAL PRODUCTS CO.,LTD

TEL 0086-574-87706193

WHATSAPP/WECHAT/TEL:008613566542343

EMAIL:SALES@HAORUNMED.COM

WWW.HAORUNMED.COM